
## Kinesiologisches Tape

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**Professional Instruction:**  
**ARTZT vitality® Kinesiological Tapes for Relief of Lumbar Spine Complaints**

### 1. Objectives

The main objective of the application of elastic tapes is to positively influence (postural) complaints. The main effect of the elastic tapes is to stimulate skin sensors. Thanks to their adhesive which is applied in a wavy pattern, the tapes adhere pleasantly to the skin and mainly stimulate local receptors. This so-called proprioceptive information (self-perception) is passed on to the brain more quickly than the "nociceptive" (harmful) information that is responsible for causing pain. Pain is thereby reduced in a natural manner.

Therefore, the tapes should be viewed as supportive self-therapy to promote the body's own healing processes. In addition to pain relief, important effects include improved blood circulation and the stimulation of lymph transportation (micro-circulation). This is why the tapes are not always effective immediately. After several hours or days, their effectiveness increases since the mechanisms described above come together. Therefore the tapes can and should be worn for several days. After approximately 5-7 days, their position should be examined and the tapes should be reapplied if necessary. More extended wear has no harmful consequences, even without examination. The tapes can be replaced as often as desired until the condition improves. Verifying that the skin is sound and unbroken is recommended prior to each application.

### 2. Target group

People with localised complaints in the spine region that are not caused by abnormal changes (e.g. herniated discs).

In case of radiating pain, especially in combination with sensitivity or motor disturbances, consulting a doctor or physiotherapist is highly recommended.

### 3. Indications – areas of application

Mild to moderate complaints in the area of the lower spine (lumbar spine).

### 4. Precautions, safety information

When the complaint is first noticed, the cause should be examined medically in order to exclude the possibility of causes that require medical treatment.

The attending doctor or a physiotherapist should be consulted in case of severe or long-term complaints.

### 5. Effect and particularities of the tapes

The main effect of the tapes is the mechanical stimulation of surface receptors. This information dampens pain perception in the central nervous system (brain) by modifying the transmission of stimuli in the spinal cord (spinal level).

The tapes can be used locally in a targeted manner at the location where the pain is perceived (symptomatic) and also "at a distance" along the pathways of the nervous system (reactive), e.g. along the spinal cord.

The tapes do not contain harmful additives and there are no known side effects (see below). The adhesive that is used has been dermatologist tested and is considered especially skin-friendly. This is why there is also normally no risk for persons allergic to conventional plasters. Combinations with other types of therapy are possible and / or recommended (supporting therapy).

Thanks to the high level of wearer comfort of the tapes, all day-to-day work and sports activities can and / or should be conducted as usual with the tapes in place. The same applies to showering. The material is breathable, permeable to moisture and liquids and dries quickly.

### 6. Application method and removal of the tapes

The skin must be dry and free of grease. Hair should be trimmed (long-hair trimmer) or gently removed entirely by shaving (if possible, the day before applying the tapes).

Intensively pressing the tapes in place and rubbing (especially along the edges) after application ensures the adhesive will adhere properly.

To remove the tapes, carefully and gradually pull in the direction of hair growth. The tapes are especially easy to remove in the showering. Skin cream or body oil help relax the skin after removing the tapes.

The lumbar spine tapes consists of a wide strip with a length of approximately 20 cm and a strip with a length of approximately 25 cm that is split in the middle (so-called Y shape). However, the length is individually adjustable. After adjusting the length, the edges have to be rounded again. Both strips are applied to the lumbar region. Because of this positioning, you need a helper in order to apply the tapes.

Take both strips out of the package and sit down on a stool or a chair without a backrest.

1) Remove the backing paper from the Y-shaped tape strip. Avoid touching the adhesive surface during this process.

2) With the end that is not split facing down towards your feet, apply the bottom end of the tape approximately two finger widths above the gluteal cleft and secure this starting point

3) In order to slightly extend the lumbar spine region, round your lower back (tilt the pelvis backwards) until you feel slight tension or reach the limit of your range of motion.

If you experience discomfort or the level of existing discomfort increases, stop the movement and straighten up until the level of discomfort decreases or the feeling goes away.

4) Now apply the split tape strips to the musculature on both sides of the spine without tension and press them into place. Make sure no folds or ripples form in the process.

5) Now carefully tear off the backing paper in the middle of the other tape strip. Avoid touching the adhesive surface during this process.

6) Stretch the centre one-third of the tape as far as possible and apply it horizontally in the area where pain is felt, centred over the "Y"

Press it into place before you apply the two outside ends of the tape without tension and press them down

Once again, make sure no folds or ripples form in the process.

### 7. Precautions, safety information

- » Always cut the tapes to size prior to application, never while you are applying them.
- » Always be careful when pulling clothing over the edges of the tapes; avoid loosening the tapes.
- » Avoid the formation of folds and ripples, especially where they could cause pressure.
- » If your complaints persist or worsen, consult your doctor and / or physiotherapist. Never apply the tapes if your complaints have not been diagnosed.
- » Remove the tapes immediately in case of intolerance.
- » Store the material at room temperature and avoid direct exposure to sunlight.

### 8. Side effects

There are no known side effects with correct application.

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH. Unauthorized use is strictly prohibited. © 2009 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

**Authors:**  
Jürgen Förster,  
Lehrphysiotherapeut  
Christian Röhrs,  
Sportphysiotherapeut DOSB  
Astrid Buscher,  
Dipl. Sportökonomin