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Unstable surfaces have been used in sports training and therapy since the 1960s. The success factors are as simple as they are appealing: ease of use is complemented by efficiency as well as a host of possible uses and areas of application. The main objective is to improve your balance. Stability trainers have recently found another useful indication in the prevention of falls. While the surface characteristics and the material of the devices are made of usually assure slip-proof contact between the body and the device, this is detrimental to the natural mobility of the human foot. In particular, the dissociation of the forefoot and hindfoot is suppressed. The Swiss medical scientist Alois Brügger pointed out the importance of "inverse torsion" for foot function back in the 1990s.

In developing a system of therapy panels, which they called Mini Stability Trainers, the objective of two French physiotherapists was to facilitate effective and efficient therapy while taking physiological aspects into account to the greatest possible extent. Recent studies confirm the superiority of the Mini Stability Trainers compared to conventional material.

Areas of application

The four therapy panels of the Mini Stability Trainers are used to improve posture (while standing) as well as movement, locomotion and / or the gait pattern. The combination of interchangeable immobile and mobile therapy panels makes it possible to complete both static and dynamic exercises. The separate use of the forefoot and hindfoot (dissociation) in a complex, natural movement pattern is unique. In addition to coordinative aspects and strength training, the main focus is on kinaesthetic perception. The differentiated perception of (joint) positions and (body) movements in space is an essential prerequisite in order to maintain and / or restore the physiologically correct axes of the legs and feet.

Harmonising the stress ratio within a muscle and between the various muscle function groups facilitates the economical use of the lower extremities on a day-to-day basis and during sports activities. This type of sensorimotor therapy lowers the actual risk of injury (prevention) and improves performance. The Mini Stability Trainers are also suitable for early functional therapy, rehabilitation and prevention.

These instructions offer suggestions and inspiration for so-called sensorimotor therapy and / or sensorimotor training (SMT) with the Mini Stability Trainers in a vivid manner with easy to understand pictures.

The (basic) exercises that are shown and the changes to the panel positions permit individual and appropriate therapy without the risk of injury.

In order to further increase performance, we have described some everyday and sports training examples that combine high efficiency with maximum safety and therapy enjoyment.

Material and characteristics

The Mini Stability Trainers are made of 100 percent washable plastic that can be disinfected. The square panels with a side length of 10 cm exhibit specific characteristics. Not only are the latter typical for the material, they also guarantee safety during therapy – with proper handling and care. The specific characteristics of the panels are also trend-setting in terms of therapy planning.

A set consists of four coloured panels. The colours serve to assign the panels to the desired sensorimotor reactions. In order to assure safety, effectiveness and efficiency, determining these specific reactions should be reserved for professionals.

A green panel with two parallel, peripheral half-rolls on the lower side is absolutely stable and does not induce any enhanced sensorimotor reactions. It is used as a fixed point and is a valuable tool to control the level of therapy and localise the sensorimotor reactions.

The yellow and blue panels, each with one central half-roll, are highly unstable and induce significant sensorimotor reactions of the targeted areas depending on

their position. They can be positioned parallel, vertically or diagonally to the foot axis. The type and intensity of the sensorimotor reactions correspond to the position of the panels. Therefore, the exercises should be completed and controlled under the guidance of an expert.

With the red panel that has a central hemisphere, the sensorimotor reaction becomes three-dimensional and correspondingly more complex.

The height of the half rolls or hemisphere was selected purposefully. It is sufficient in order to trigger therapeutically beneficial and meaningful reactions of the sensorimotor system. On the other hand, the equilibrium organ is not significantly stimulated by the range of motion and the subsequent minor change of the body position in space. Through this purposeful focus, the therapeutic benefits are maximised while the risk of injury is minimised.

Handling

Ensure that the Mini Stability Trainers are used on a flat, firm surface. Using the panels on slipper surfaces should be avoided for safety reasons.

Please note that the tilting angle and range of motion of the mini stability trainers are affected by softer floor.

The position of the panels depends on the therapy objective and is therefore decisive in terms of therapy success and safety.

The diagnosis made by a physician and / or the findings of a physiotherapist form the basis for the efficient and safe use of the Mini Stability Trainers in therapy.

General counter-indications, such as signs of local inflammation, also need to be taken into account in the course of therapy with the Mini Stability Trainer panels.

- » In principle, all exercises must be conducted so that they are pain-free.
- » All movements must be controlled at all times; they are completed slowly at first and more quickly later on.
- » Full breaks should be taken between repeats and series in order to assure high-quality therapy.
- » Quality is more important than quantity.

Initial position

The initial position has a key influence on the therapy results (outcome). In principle, the initial position is as active as possible and as passive as necessary. The general body posture should be as upright and economical as possible at that moment. The correction of posture by the therapist prior to therapy may be advisable.

While the basic exercises are completed in a standing position or while standing on one leg, other initial positions are also possible. The choice depends on the therapy objectives and control. Regardless of the body position in space, a correct posture should be maintained in the initial position.

Duration

The static exercises (poses) are held for 20 to 30 seconds at a time in 1-3 series.

Repetitions

The dynamic exercises (movements) are repeated 10 to 15 times in 1 to 3 series.

Progression

In addition to the usual stress parameters (intensity, scope, frequency and duration), the following progressions apply:

- » From simple to complex.
- » From slow to fast.
- » From active (conscious) to reactive (subconscious).
- » With visual control, then without.

The stress parameters listed below are relevant at the functional level (perception, control, cybernetics). Therapy is continued or intensified in case of positive parameters and reduced or discontinued in case of negative parameters.

Parameter	Positive	Negative
Strength	Increase	Decrease
Range of motion	Increase	Decrease
Coordination	Improvement	Deterioration
Pain	Decrease	Occurrence or increase
Rigor	Decrease	Occurrence or increase
Subjective fatigue	None	Occurrence

Training information

The four Mini Stability Trainer panels facilitate a variety of exercises. In addition to the basic exercises, special training variations with various degrees of difficulty and different initial positions are also described. The position of the feet is extremely important and must be assumed before each exercise; it should be verified continuously.

The natural (leg) foot axes serve as the initial position:

- » The foot is flat on the floor.
- » The long toe muscles, i.e. the extensors, are relaxed according to the load.
- » The arches of the foot are defined in a physiological manner.
- » The weight on the legs is evenly distributed between the forefoot, hindfoot and arch of the foot.

Most exercises are completed under (full) load, i.e. in a standing position. Nevertheless, it may be expedient to begin therapy from a sitting position or while kneeling on one knee in the beginning or when training under partial load.

Exercises

Basic exercises

The basic training series consists of 21 select exercises and forms the introduction to sensorimotor therapy with the Mini Stability Trainer panels. They are arranged progressively according to the level of difficulty and demand.

The series consists of exercises that respectively train the hindfoot or forefoot in particular and exercises that mainly train the hindfoot or the forefoot as well as exercises that train both the hindfoot and the forefoot equally. Differentiated strain on specific muscle groups rounds out the training series.

The following table presents a summary of the basic training series:

	Hindfoot	Forefoot	Differentiated Strength Training	Number of Exercises
Basic exercises I	+	-	-	4
Basic exercises II	-	+	-	4
Basic exercises III	+	+/-	-	4
Basic exercises IV	+/-	+	-	4
Basic exercises V	+	+	-	3
Basic exercises VI	-	-	+	2

Basic exercises I:

specific hindfoot training

Basic exercises II:

specific forefoot training

Basic exercises III:

predominant hindfoot training (with moderate forefoot training)

Basic exercises IV:

predominant forefoot training (with moderate hindfoot training)

Basic exercises V:

balanced hindfoot and forefoot training

Basic exercises VI:

differentiated training of specific muscle function groups (eversors, invertors)

All basic exercises may be varied by changing the arm position (on the sides or crossed) and / or visual control (eyes open or closed).

Mastering one level is essential before moving on to the next one! Increasingly shifting the weight of the body and / or leg to the affected leg (from approximately 50 to 100 percent) represents another progression option.

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Position the yellow panel diagonally under the hindfoot (at an angle of 45 degrees to the longitudinal axis of the foot) and the green panel under the forefoot. The toes should not project beyond the edge of the panel. Both panels are centred.

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Position the yellow panel diagonally under the hindfoot (at an angle of 45 degrees to the longitudinal axis of the foot, but reversed compared to exercise 1) and the green panel under the forefoot. Both panels are centred.

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Position the yellow panel under the hindfoot so the half-rolls are parallel to the longitudinal axis of the foot, and the green panel under the forefoot. Neither the toes nor the heel should project beyond the edge of the panel. Both panels are centred.

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Position the red panel so it is centred under the hindfoot and the green panel under the forefoot. Neither the toes nor the heel should project beyond the edge of the panel.

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Position the blue panel diagonally under the forefoot (at an angle of 45 degrees to the longitudinal axis of the foot) and the green panel under the hindfoot. The heel should not project beyond the edge of the panel. Both panels are centred.

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Position the blue panel diagonally under the forefoot (at an angle of 45 degrees to the longitudinal axis of the foot, but reversed compared to exercise 5) and the green panel under the hindfoot. Both panels are centred.

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Position the blue panel under the forefoot so the half-rolls are parallel to the longitudinal axis of the foot, and the green panel under the hindfoot. Neither the toes nor the heel should project beyond the edge of the panel. Both panels are centred.

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Position the red panel so it is centred under the forefoot and the green panel under the hindfoot. Neither the toes nor the heel should project beyond the edge of the panel.

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**Predominant hindfoot training
(with moderate forefoot training)**

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Position the blue panel diagonally under the forefoot (at an angle of 45 degrees to the longitudinal axis of the foot) and the yellow panel under the hindfoot so that the half-rolls are parallel to the longitudinal axis of the foot. The heel should not project beyond the edge of the panel. Both panels are centred.

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Position the blue panel diagonally under the forefoot (at an angle of 45 degrees to the longitudinal axis of the foot, but reversed compared to exercise 9) and the yellow panel under the hindfoot so that the half-rolls are parallel to the longitudinal axis of the foot. Both panels are centred.

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**Predominant forefoot training
(with moderate hindfoot training)**

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Position the yellow panel diagonally under the forefoot (at an angle of 45 degrees to the longitudinal axis of the foot) and the blue panel under the forefoot so that the half-rolls are parallel to the longitudinal axis of the foot. Neither the toes nor the heel should project beyond the edge of the panel. Both panels are centred.

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Position the yellow panel diagonally under the hindfoot (at an angle of 45 degrees to the longitudinal axis of the foot, but reversed compared to exercise 13) and the blue panel under the forefoot so that the half-rolls are parallel to the longitudinal axis of the foot. Both panels are centred.

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Position the yellow panel diagonally under the hindfoot (at an angle of 45 degrees to the longitudinal axis of the foot) and the red panel under the forefoot. The toes should not project beyond the edge of the panel. Both panels are centred.

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Position the yellow panel diagonally under the hindfoot (at an angle of 45 degrees to the longitudinal axis of the foot, but reversed compared to exercise 11) and the red panel under the forefoot. The toes should not project beyond the edge of the panel. Both panels are centred.

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Position the blue panel under the forefoot so that the half-rolls are parallel to the longitudinal axis of the foot and the yellow panel under the hindfoot so that the half-rolls are parallel to the longitudinal axis of the foot. Neither the heel nor the toes should project beyond the edge of the panel. Both panels are centred.

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Position the red panel under the forefoot and the yellow panel under the hindfoot so that the half-rolls are parallel to the longitudinal axis of the foot. Neither the heel nor the toes should project beyond the edge of the panel. Both panels are centred.

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Position the red panel under the hindfoot and the blue panel under the forefoot so the half-rolls are parallel to the longitudinal axis of the foot. Neither the heel nor the toes should project beyond the edge of the panel. Both panels are centred.

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Differentiated training of the eversors

Position the blue panel under the forefoot so the half-rolls are parallel to the longitudinal axis of the foot, but slightly offset towards the inside, and the green panel under the hindfoot. Neither the toes nor the heel should project beyond the edge of the panel. The green panel is centred.

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Differentiated training of the invertors

Position the blue panel under the forefoot so the half-rolls are parallel to the longitudinal axis of the foot, but slightly offset towards the outside, and the green panel under the hindfoot. Neither the toes nor the heel should project beyond the edge of the panel. The green panel is centred.

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The special exercises offer a progression of the basic exercises. They follow the principles of sensomotoric training and turn out to be an additional improved performance.

The exercises are performed in 1 to 3 series, being held statically 20 to 30 seconds, or repeated dynamically 10 to 15 times.

As already mentioned the initial position and the body posture are essential elements for the success of the therapy. The upright posture should be controlled before each exercise. The panels must be positioned as described under the basic exercises and adapted to the objectives and the abilities of the patient.

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Maintain an upright standing position.

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Stand on one leg, free leg with floor contact (toes).

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Stand on one leg, free leg without floor contact (behind the body).

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Stand on one leg, free leg without floor contact (in front of the body).

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Stand on one leg, free leg with floor contact, arms to the sides; practice with eyes open first, subsequently with eyes closed.

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Stand on one leg, free leg with floor contact, arms crossed in front of the chest; practice with eyes open first, subsequently with eyes closed

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A flat, slip-proof surface is particularly important for the following dynamic exercises!

The series of exercises presented here represents the way the foot absorbs stress in the first half of the gait cycle: First heel impact, then the progressive shift of weight as the heel of the rear leg comes off the floor.

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Everyday movements when standing up and sitting down are made more effective through the use of the Mini Stability Trainers. The selection of Mini Stability Trainers must be individually established. It has to match the capabilities of the patient and the objectives. Movement transitions are conducted with the feet parallel in the walking position.

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By positioning the Mini Stability Trainers under one or both feet, the straddle position is made truly challenging for the sensorimotor system. In addition, the weight of the body can be shifted between the feet!

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
From a secure standing position, the body weight is first transferred to the forefoot and then the ankle is extended. The other foot has floor contact and keeps the position steady!

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The effectiveness of the familiar “functional squat” can be significantly increased with the use of the Mini Stability Trainers.

UPGRADE YOUR VITALITY!

ARTZT
vitality

Produktprogramm ARTZT vitality® | Product range ARTZT vitality® Medical Fitness - Made in Germany

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- Flowin® Personal Training
- Gymstick™ Health
- BOSU® Balance Training

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ARTZT vitality® is a registered trademark of Ludwig Artzt GmbH. Unauthorized use is strictly prohibited.
© 2009 Ludwig Artzt GmbH. All rights reserved. Printed in Germany.

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der Ludwig Artzt GmbH. Unerlaubter Gebrauch wird strafrechtlich verfolgt.
© 2009 Ludwig Artzt GmbH. Alle Rechte geschützt. Gedruckt in Deutschland.